# Z-12 Knee Brace

Application Instructions CP020179 Rev A 12/08

## 1. Aprire le fibbie

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## 1. Unfasten Buckles

Release all four straps on the back side of the brace by squeezing the levers on either side of the buckle and withdrawing the strap tab. Loosen the front strap(s).

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## 3

## 3. Position the Brace

With the brace still in the prefit position, locate the top of the knee cap. Position the center of the hinge on an imaginary line that goes across the top of the knee cap. Push the hinges posteriorly until the hinges are in the center of the knee joint when viewed from the side.

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## 4

## 4. Anchor the Brace

Wrap the upper rear calf strap around the calf as high as possible over the calf muscle.

Insert the over-center tab into the base. Push the tab until you hear the hinge click on both sides. Adjust the strap by releasing the tab and pull the strap until a snug and comfortable fit is achieved.

All straps have additional length adjustment at the opposite end of each strap. See step 6.

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## 6

## 6. Adjust Strap Length

All straps have additional length adjustment at the end opposite the tab. Release the strap and adjust the length of the strap using the hook and pile fastener.

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**7. Fasten Rear Thigh Straps**  
 Following the same instructions as in Step 5, fasten the lower and then upper rear thigh strap.

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
**8. Adjust Front Tibia Strap**  
 Tighten the front tibia strap. If the strap overlaps too far, see step 6 for adjustment instructions.

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
Centre la almohadilla de la correa.

**9. Adjust Front Thigh Strap (PCL/Combined Instability Versions)**  
 Tighten the front thigh strap. Note: Overtightening the front thigh strap can cause the brace not to function properly.

Stand up and walk or jog around to see how your new AXIOM or AXIOM-D brace feels. Straps should not be so tight that blood circulation is restricted. Re-adjust as needed. If this is your first time to wear a functional knee brace, it might take a few hours of wear to get used to how it feels.

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## Patient Information

### Special Adjustment for Anterior Cruciate Ligament Injuries

In order to provide increased resistance against tibial subluxation (forward movement of the lower leg relative to the upper leg) when the leg is extended or to reduce force on a reconstructed ACL, it may be desirable to pretension the straps to provide a "shear force" across the knee by increasing the tension on the front calf strap on the tibia and the lower rear thigh strap on the back of the thigh. All straps for ACL patients should be adjusted with the leg bent at 60°-80°. Loosen the front thigh strap and increase the tension on the lower rear thigh strap and, if necessary, the front shin strap to provide a pre-loaded strap force that helps prevent anterior tibial subluxation.

### Special Adjustment for Posterior Cruciate Ligament Injuries

In order to provide increased resistance against posterior tibial subluxation (backward movement of the lower leg relative to the upper leg) as the leg bends or to reduce force on a PCL reconstruction, it may be desirable to pretension the straps to provide a "shear force" across the knee by increasing the tension on the upper rear calf strap and the front thigh strap. All straps for PCL patients should be adjusted with the leg in full extension. Loosen the lower rear thigh strap and the front calf strap and increase the tension on the front thigh strap and the upper rear calf strap to provide a pre-loaded shear force that prevents posterior tibial s ubluxation.

NOTE: If needed, apply the included spacers behind the condyle pads to achieve proper condular pressure on the medial and/or lateral side of the knee.

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This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient. This device is generally not intended to prevent injuries or reduce or eliminate the risk of reinjury.

**Indications:** For increased support for patients who have continued symptoms of giving way, poor quadriceps or hamstring strength, or a desire for an early return to sports activities. Indicated for increased support of knees following injury to or reconstruction of the anterior or posterior cruciate ligament, medial or lateral collateral ligament, or menisci either alone or in combinations when properly adjusted.

**Contraindications:** Contraindicated in patients with triple varus knees that have damage to the postero-lateral structures combined with stretched lateral structures and varus thrust.

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**WARNING:** Do not leave this device in the trunk of a car on a hot day or in any other place where the temperature may exceed 150° F (65° C) for any extended period of time. **If any additional pain or symptoms occur while using this device, discontinue use and seek medical attention.**

**Warranty:** This device is warranted to be free from defects in material and workmanship for a period of 90 days for wraps, pads and straps and normal wear components; 12 months on all other parts; and 3 years against hinge breakage. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect. This device is expressly subject to Bledsoe's Product Line Limited Warranty. Copies of this limited warranty can be obtained from your medical professional, via the internet at www.BledsoeBrace.com, or by calling 1.800.527.3666.

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