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www.djortho.com
Fax: (800) 936-6569
Toll Free: (800) 936-6569
Visita, CA 92083-8339 U.S.A.
2985 Scott Street
dj Orthopedics, LLC

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DEVELOPED WITH DR. TOM SAWA
DONJOY SHOULDER STABILIZER

dj Ortho and DonJoy are trademarks of dj Orthopedics, LLC
dj Ortho y DonJoy es una marca registrada de dj Orthopedics, LLC
dj Ortho et DonJoy est une marque déposée de dj Orthopedics, LLC
dj Ortho und DonJoy ist ein eingetragenes Warenzeichen von dj Orthopedics, LLC
dj Ortho e DonJoy sono marchi depositati dalla dj Orthopedics, LLC

13-3110-0-0000

© dj Orthopedics, LLC


APPLICATION INSTRUCTIONS

Before applying the brace, please read these instructions completely and carefully.
Correct application is vital to the proper functioning of the brace.


INDICATIONS: • Help reduce and minimize further injury to an already existing shoulder injury.

BRACE APPLICATION

A Apply the brace over the injured shoulder. Grasp one end of brace at bottom and pull over hook and loop closure, making sure to close bottom tighter than top. (Figure 1) The fit should be snug but should not restrict breathing.


B Secure the hook and loop strap through the "D" ring to lock in side closure. Attach the hook and loop closures at the biceps tying to get as tight a fit as possible, this will minimize any rotation within the biceps cuff. (Figure 2)


ABDUCTION CONTROL

C To achieve simple abduction control as in figure 3, after applying the brace, attach the glenohumeral straps to the biceps cuff with the stitching (see figure 6) attaching the front and back straps positioned as in figure 4a. This positioning is very important.


IMPORTANT

The positioning of the stitching is critical to maximizing the allowable forward flexion while at the same time minimizing external rotation. In figure 4 as you move the positioning from A to C and shorten the front strap you will minimize external rotation.


D Simply tighten or loosen the straps through the "D" rings and hook and loop to achieve the desired range of motion and abduction control. (Figure 5)


ADVANCED FITTING

E The shoulder brace is very effective in maximizing abduction in forward flexion. (Figure 7) This is ideal for people who must reach over their head but still require prevention of external rotation and backward extension in abduction.


F Loosen the back strap, if necessary take the back strap completely out of the "D" ring. Position stitching in front of biceps as in figure 4. It is now important to tighten the front strap to control backward extension. (Figure 8)

G A/C Separations

The shoulder brace is also effective in supporting acromioclavicular separations. If the joint is still tender, a donut or protective pad can be placed under the brace. Tighten the A/C strap by adjusting the strap through the front and back "D" rings. The long or medium auxiliary straps can also be threaded through the loops over the A/C and attached to the hook and loop at the front and back of the brace. (Figure 9)

USE AND CARE

- Machine wash, cold water. AIR DRY.

WARNING: IF YOU EXPERIENCE ANY PAIN, SWELLING, SENSATION CHANGES, OR ANY UNUSUAL REACTIONS WHILE USING THIS PRODUCT, CONSULT YOUR MEDICAL PROFESSIONAL IMMEDIATELY.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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